

e-PAGEL Instruction Manual

July 20th, 2016 Ver.2

1. Safety precautions

Before using the product, read this manual thoroughly at first. Do not start the operation until you understand the contents of manual. This document explains only methods utilized for specified purposes. Do not use the product for any purpose or by any method not described in this manual. If it is used for any purpose or by any method not described in this manual, an operator should take responsibility for all required safety measures and contingencies. Also, read a manual of equipment used with it simultaneously.

2. Application purpose

「e-PAGEL」is a precast polyacrylamide gel for electrophoresis of protein and DNA. An electrophoresis device (precaster gel specification) exclusively for ATTO Mini-gel size are required.

3. Package

Product name	Size	Package
e-PAGEL	Gel size	10/pk
	90(W)x83(H)x1mm(t)	
	Glass plate size	
	120(W)x100(H)x2mm(t)(total 5mm)	

4. Components

Product name	Major components
e-PAGEL	Polyacrylamide gel

This product doesn't include a notifiable material exceeding to regulated amount for exclusion decided by PRTR Law, Poisonous and Deleterious Substances Control Act, and Industrial Safety and Health Law.

5. Preservation method

- Please keep refrigerated. Quality is impaired when freezing.
- The expiration date is displayed on the outer box and gel packaging bag.

6. Disposal method

- Follow a disposal method decided by the organization you belong to.
- Materials Plate: Glass Packaging: PET nylon

7. Necessary things other than the product

- Electrophoresis apparatus for ATTO Mini-gel size
- Power supply (recommended output over 300V, 150mA)
- Electrophoresis buffer etc.

8. Precautions for use

- Please keep refrigerated. (5~10°C) Please do not place near the cold air outlet. It may freeze.
- Freezing causes deformation such as bubbles, plate peeling, swelling and shrinkage, and can not be used.
- Please keep upright according to the mark on the outer box. If not, the quality may deteriorate.
- Follow the mark and store the gel with the well facing down. Please note that storing in the wrong direction will cause the quality to deteriorate.
- Please open the bag just before use. After opening, please use it immediately as the quality deteriorates.

- Handling the gel with bare hands may cause injury. Please wear gloves and protective clothing when handling.


9. Usage

9-1. Preparation of gel and electrophoresis buffer

1. Open the packaging bag and take out the gel.


※When it is difficult to cut, please cut it with scissors.

※If you pull out the gel forcibly, the glass will peel off from the gel. Please avoid it.



2. Pull out the comb slowly.

※Place your fingers on the convex parts of the comb and remove the comb slowly by aligning the left and right sides. Remove the comb without bending or cutting the wells.



3. Prepare the electrophoresis buffer corresponding to the sample treatment (with or without SDS treatment) and the gel type.

SDS treated protein : EzRun (Tris/Gly/SDS)

: EzRun MOPS (Tris/MOPS/SDS)

SDS-free protein and nucleic acid

: EzRun TG (Tris/Gly)

Tricine gel for low molecular weight **P-T16.5S/P-R16.5S**

: EzRun T (Tris/Tricine/SDS)

4. Wash wells with electrophoresis buffer.

9-2. Electrophoresis

1. Place the gel in an electrophoresis device exclusively for ATTO Mini-gel size and pour the electrophoresis buffer.

※Please set the gel according to the instruction manual attached to the electrophoresis device.

2. Apply the appropriate amount of sample to each well.

※The maximum applied amount is expressed as about 60% of the maximum capacity of the well.

Comb (Number of wells)	Well size	Maximum applied amount
14 wells	4.2(W)×10(H) mm	24 μL
18 wells	2.9(W)×10(H) mm	18 μL
2-D	78(W)×25(H) mm	AgarGEL 1ea

3. Set the power supply referring to the following table.

Electrophoresis buffer		Voltage	Current	Time
EzRun Tris/Glycine/SDS	C.C.	Starting: 75-90V Ending: 180-220V	20mA/gel C.C.	75-80min
	C.V.	250V C.V.	Starting: 70-100mA Ending: 30-50mA	25-35 min
EzRun MOPS Tris/MOPS/SDS	C.C.	Starting: 65-85V Ending: 130-180V	20mA/gel C.C.	60-70min
	C.V.	Starting: 100-120V Ending: 230-250V	20mA/gel C.C.	60-90 min
EzRun TG (Nucleic acid) Tris/Glycine	C.C.	Starting: 65-75V Ending: 130-160V	20mA/gel C.C.	130-150min

※C.C.; Constant Current, C.V.; Constant Voltage

※In case of C.V./C.C., set the allowable voltage and current value refer to the table.

※The current or voltage value (measured value) other than the setting values are the approximate. It depends on the gel concentration.

※For constant voltage setting, please set it to 250V regardless of the number of gels.

For constant current setting, set the current value as “the number of gels × 20 mA”.

4. Apply the power and perform the electrophoresis.


※After electrophoresis, buffer solution can be hot. In case of 20 mA/gel constant current setting, it may increase to around 25°C, and in case of 250 V constant voltage setting, it may increase to 35-50°C.

9-3. End of electrophoresis

1. Stop migration when the dye front reaches the 5~10 mm upper position from the lower end of the gel.

※Please stop the operation according to the instruction manual attached to the electrophoresis device.


2. Remove the gel plate from the electrophoresis tank, insert a flat object such as a spatula between the glass, and gently move it up and down to open the gel plate. Remove the top glass plate.



3. Cut between the gel and the spacer with a spatula or scalpel.

※Moist spatula moves smoothly and prevents damage to the gel.

4. Transfer the gel to the vat containing the staining solution. Hold the gel plate with the gel side down, insert a spatula between the gel and the gel plate, and peel off the gel.



5. Slightly shake the vat to immerse the whole gel.



※CBB/silver/reverse/fluorescent staining, etc. can be use.



※When staining **P-T/R16.5S** with **EzStain AQUA**, use a microwave to prevent the elution of low molecular weight proteins. If you using a hand-made CBB staining solution, soak in the 50% methanol or 10% acetic acid for 20~30 min before staining with the CBB to fix the proteins.


9-4. 2-D electrophoresis

1. Dip molecular weight marker (~3 µL) on 3 mm×3 mm filter paper. Place it on the edge of the gel in contact with the upper surface of the 2-D gel.
2. Place the equilibrated 1-D gel in contact with the upper surface of the 2-D gel.



※Make sure that the 1-D and 2-D gels are in close contact. A gap between the gels may cause the electrophoresis to be disturbed.

3. Drop 1% agarose solution and attach the 1-D and 2-D gels.



※Please use the agarose solution after cooling to around 60°C

4. Perform electrophoresis (see section 9-2).

※Please follow the instruction manual attached to the agarGEL and electrophoresis device.

10. Reference

10-1. Gel concentration and separation range

Electrophoresis buffer

Protein:

EzRun (Tris/Gly/SDS)

Nucleic acid:

EzRun TG (Tris/Gly)

P-T/R16.5:


EzRun T (Tris/Tricine/SDS)

EzRun MOPS:

See section 10-2

Gel conc.	M.W. range (Protein)	M.W. range (Nucleic acid)
7.5%	40-400 kDa	200-3000 bp
10%	25-300 kDa	100-2000 bp
12.5%	10-250 kDa	70-1800 bp
15%	5-200 kDa	50-1500 bp
5-20%	5-400 kDa	30-2500 bp
10-20%	5-300 kDa	30-2000 bp
E-T/R16.5S	1-75 kDa	-

10-2. Electrophoresis pattern



0,000 001 = 10⁻⁶ micro
0,000 000 001 = 10⁻⁹ nano
0,000 000 000 001 = 10⁻¹² pico
0,000 000 000 000 001 = 10⁻¹⁵ femto
0,000 000 000 000 000 001 = 10⁻¹⁸ atto



(주)아토코리아

Providing research solutions for
Biochemistry / Molecular Biology /
Genetic engineering development

대전광역시 유성구 복용로40번길 5-6, 102호 (복용동)

■ Tel. (042) 822-1117

■ URL www.attokorea.co.kr

■ Fax. (042) 825-1116

■ E-Mail info@attokorea.co.kr